


How to fit your stick clip (HX15153R_1)

Keep your reacher/grabber with you wherever you go.




This stick clip is specifically designed to fit snugly on all 'Handi-Grip' & 'ArthriGrip' reachers

To fit on to your reacher:



***Pinch the 'rabbit ears'
to open the clip***



***Push onto the side
of the reacher bar***



***Make sure it 'clunk-
clicks' into place***

*Clip it on to a crutch, chair leg
walking pole/stick/frame or
your wheelchair – simple.*

*For more information
and other useful products
please visit:
www.HHADLEssentials.co.uk*



The Helping Hand Company
Bromyard Road, Ledbury
Herefordshire HR8 1NS
Tel: 01531 635678
Fax: 01531 635670
shop@helpinghand.co.uk


**The Helping Hand
Company**